

FCC Markings

WARNING: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

SUPPLIER'S DECLARATION OF CONFORMITY

Brand: ConairCare®

Model No.: CN600

Description: Bathroom Scale

Responsible Party: Conair LLC

Standards: FCC Part 15 Subpart B, ANSI C63.4-2014

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) This device must accept any interference received, including interference that may cause undesired operation.

DATE OF ISSUE: July 04, 2023

Please register this product at register.conair.com


FSA/HSA

In order to receive reimbursement, you may have to provide a Letter of Medical Necessity (available by scanning the QR code here) completed and signed by your doctor.



Register Your Product Now

Register in two easy steps to receive premium customer service, special offers, and more.



Step 1

Take a photo of the camera icon at left

Step 2

Send it in using your mobile phone

Smartphone (US ONLY): text photo to 71403

You can also register via the web: register.conair.com

Need Help? Visit photoregister.com/help or text HELP to 71403.

Important: Failure to complete and submit this form does not diminish your warranty rights.

Escanear para español

Escanear este código para ver la versión en español del embalaje y del manual de instrucciones.



Service Center
Conair LLC
7475 N. Glen Harbor Blvd.
Glendale, AZ 85307

Questions
about your scale?
800-326-6247
www.conair.com



©2023 CONAIR LLC, GLENDALE, AZ 85307
MADE IN CHINA

23AD088744 **IB-18158**

HOME HEALTHCARE HIGH-CAPACITY SCALE



CONAIRCARE®
PRECISE. SAFE. FAST.

Instructions for use of the CN600

*To get the greatest benefit from your purchase,
please read this entire brochure before using your scale.*

PRODUCTION INKS



Black

All prints 1 color: Black

Other instructions:

Hi res in place

Bleed edges


INSTRUCTIONS FOR USE

How to Operate

1. Remove all packing materials.
2. Remove the isolator tab from the battery compartment – without opening the battery cover at the bottom of the scale, gently pull out the isolator tab. (If the tab cannot be pulled out or if the display does not turn on after the tab is pulled, remove the cover. Gently pull out the isolator tab. Ensure the batteries are installed correctly. Replace the battery cover.)
3. For best accuracy, and consistent readings, place scale on tile or hardwood floor, rather than uneven, flexible, or soft flooring (such as carpet).
4. Stand on platform without shifting or moving until weight is displayed (4–5 seconds). **NOTE:** Digital display automatically turns off after 20 seconds.

Selecting Units of Measurement

When scale display is on, press KG or LB button on the underside of the scale for measuring in kilograms (kg) or pounds (lb).




Talking Mode

The optional talking mode clearly announces weight readings in either English or Spanish for those who are visually impaired.

When scale display is on, press the speaker button on the underside of the scale to select a talking mode: English (EN), Spanish (ES) or Silent (X).

Caregiver Display Mode

The patent-pending Caregiver Display Mode allows a caregiver to easily read the display while facing the person being weighed.



When scale display is on, press the button above the LCD display to activate the caregiver flip display feature. Press again to return to normal display.

Helpful Hints

Your ConairCare® scale is a precise measuring instrument that is most accurate when weighing a stationary object. To ensure the most accurate readings, always try to stand on the same area of the scale platform and **DO NOT MOVE**. Should the weight on the scale exceed the scale's capacity (660 lb.), you may see an error message — "Err."

Electronic sensors are sensitive. Be careful not to drop or jar the scale. Place it gently onto floor surface, and store where it will be protected from impact. The scale is an electronic instrument and should never be submerged in water. Clean with a damp cloth only.

"Lo" will appear on the display if the batteries are low and need to be replaced. Remove and replace batteries. Clean the battery contacts and also those of the device prior to battery installation. Ensure the batteries are installed correctly with regard to polarity, then replace the battery cover again. Remove used batteries promptly. When replacing the batteries, make sure you dispose of the old batteries in accordance with applicable regulations. They should be disposed of in a location specifically designated for that purpose so they can be recycled safely and without posing any risk to the environment. Do not burn or bury them. See Battery Notes and Battery Safety Information for more important info on batteries.

Do not attempt to open the scale or to remove any components. Servicing should be done by qualified technicians only. See warranty information below.

LIMITED TEN-YEAR WARRANTY (US and Canada only)

Conair will repair or replace, at our option, your appliance free of charge for 120 months from the date of purchase if the appliance is defective in workmanship or materials.

To obtain service under this warranty, return the defective product to the service center listed, together with your sales receipt. California residents need only provide proof of purchase and should call 1-800-3-CONAIR for shipping instructions. In the absence of a receipt, the warranty period shall be 120 months from the date of manufacture.

ANY IMPLIED WARRANTIES, OBLIGATIONS, OR LIABILITIES, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE, SHALL BE LIMITED IN DURATION TO THE 120-MONTH DURATION OF THIS WRITTEN, LIMITED WARRANTY. Some states do not allow limitations on how long an implied warranty lasts, so the above limitations may not apply to you.

IN NO EVENT SHALL CONAIR BE LIABLE FOR ANY SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES FOR BREACH OF THIS OR ANY OTHER WARRANTY, EXPRESS OR IMPLIED, WHATSOEVER. Some states do not allow the exclusion or limitation of special, incidental, or consequential damages, so the above limitation may not apply to you. This warranty gives you specific legal rights, and you may also have other rights, which vary from state to state.

THINGS YOU SHOULD KNOW ABOUT MANAGING YOUR WEIGHT

Diet, exercise, and determination all play a part in effective weight management. Your scale is the tool that lets you chart your progress and measure your success, whether your goal is weight loss, weight gain, or maintaining your current weight.

Because weight varies, as explained below, focus on the progression and sustained, long-term changes rather than on exact, day-to-day readings. When you're dieting and exercising, small, fluctuating weight gains and losses are common.

Common Causes of Weight Fluctuations (during the day and day-to-day)

- Water retention and weight gain from consumption of large meals.
- Water retention from salt and carbohydrate intake.
- Muscle buildup (muscle weighs more than fat) from strenuous resistance training.
- Weight loss caused by dehydration due to illness, low fluid intake, or strenuous exercise.
- Temporary weight gains and losses from menstrual cycle.

Weight Maintenance vs. Weight Loss/Gain

Get into the habit of weighing yourself no more than once a week when you're on a weight management program. Weighing in at consistent intervals will give you the most accurate picture of your progress, eliminating many short-term, unrelated weight changes.

Follow These Simple Rules:

- Weigh yourself once a week.
- Weigh yourself on the same day of the week.
- Weigh yourself at the same time of day.
- Weigh yourself wearing the same clothes.
- Weigh yourself in the same place.
- Weigh yourself on the same scale.

Weight fluctuates more when a person is overweight or underweight. When the goal is to maintain current weight, weighing yourself more than once a week is a good way to stay on track.

Remember: A scale is an essential tool in any weight management program, but it only records your progress. Exercise, diet, and determination are what you need to reach your goal.

Battery Notes:

- Always purchase the correct size and grade of battery most suitable for the intended use.
- Replace all batteries of a set at the same time.
- Do not mix old and new batteries. Do not mix alkaline and standard (carbon-zinc) batteries. Do not use rechargeable (Ni-Cad, Ni-MH, etc.) batteries.
- Clean the battery contacts and also those of the device prior to battery installation.
- Ensure the batteries are installed correctly with regard to polarity (+ and -).
- Remove batteries from a scale that is not to be used for an extended period of time.

Dispose of the old batteries properly. Remove batteries from equipment that is not to be used for an extended period of time.

When batteries of different brand or type are used together, or new and old batteries are used together, some batteries may be over-discharged due to a difference of voltage or capacity. This can result in venting, leakage, and rupture and may cause personal injury.

Do not attempt to open the scale or to remove any components. Servicing should be done by qualified technicians only. See warranty information in this brochure.

Battery Safety Information

1. This product contains AA alkaline batteries.
2. Since batteries are sensitive to high temperatures, when storing them, keep in a cool, dry place.
3. Battery contacts need to be clean, both on the battery and in the appliance, before installation.
4. Batteries should be removed from the appliance when not being used for long periods of time (such as months).
5. Never recharge non-rechargeable batteries.
6. Be careful and install (+) and (-) on batteries to the correct plus and minus on device.
7. Replace all old batteries at one time. Do not mix old and new.
8. Do not mix alkaline batteries with standard or rechargeable batteries.
9. Remove discharged batteries immediately to avoid leakage.
10. Never put batteries into a fire.
11. Keep batteries away from children, especially the batteries that are small enough to ingest.
12. If a battery is ingested, seek immediate medical attention. Call the National Battery Ingestion Hotline at (800) 498-8666.
13. Never open, crush, puncture or disassemble batteries.
14. Dispose of batteries according to any state and local regulations.

Recyclable Batteries

Alkaline batteries contain various materials that must be recycled and disposed of properly. Do not dispose of in the municipal waste system. Under various state and local laws, it is illegal to dispose of these batteries in the municipal waste stream. The EPA-certified RBRC™ battery recycling seal on the alkaline batteries indicates Conair is voluntarily participating in an industry program to collect and recycle these batteries at the end of their useful life, when taken out of service in the United States or Canada. The RBRC™ program provides a convenient alternative to placing the used alkaline batteries in the trash or the municipal waste stream, which may be illegal in your area. Please call 1-800-8-BATTERY™ for information on alkaline batteries. Check battery recycling requirements and disposal bans/restrictions in your area. Conair's involvement in this program is part of our commitment to preserving our environment and conserving our natural resources. Please call 1-800-8-BATTERY™ or refer to the Call2Recycle link below for the nearest drop-off recycling location: www.call2recycle.org/locator/