
# CONAIR® Flex

## Pain Relief

## User Guide


**INCLUDES IMPORTANT SAFETY INFORMATION. READ ALL INSTRUCTIONS BEFORE USING DEVICE. SAVE THESE INSTRUCTIONS.**

Thank you for your purchase. We want you to be completely satisfied with the quality and performance of this product. Our goal is to discover new products for your health. We strive to provide exceptional customer service.

Conair Corporation

### I. Product Description

The Conair® Flex Pain Relief is a flexible, infrared light device for therapeutic healing. It delivers infrared light to tissues resulting in therapeutic effects. The device contains infrared 880nm and Red 660nm LEDs. Treatment time is controlled by the operator. Recommended treatment time should last 14 - 25 minutes on any one particular part of the body.

### II. The Conair® Flex Pain Relief Indications for Use

The Conair® Flex Pain Relief is intended for the relaxation of muscles and relief of muscle spasms; relief of minor muscle and joint aches, pain and stiffness associated with arthritis, and to increase local blood circulation.

## IMPORTANT SAFETY INSTRUCTIONS DO NOT DESTROY READ ALL INSTRUCTIONS BEFORE USING

**DANGER** – To reduce the risk of burns, electric shock and fire, this product must be used in accordance with the following instructions:

1. READ ALL INSTRUCTIONS BEFORE USING.
2. DO NOT USE WHILE SLEEPING. KEEP OUT OF REACH OF CHILDREN; DO NOT USE ON CHILDREN.
3. The device is not to be used on or by an invalid, a paralyzed person, a sleeping or unconscious person, a person with diabetes, or a person with poor blood circulation. – Do not use on areas of sensitive skin and areas that do not have sensation.
4. Temperatures sufficiently high to cause burns may result from improper use of this device. Check skin under device frequently to avoid burning and blistering.
5. Do not use in an oxygen-enriched environment or near equipment that stores or emits oxygen.
6. Do not use pins or other metallic means to fasten device in place under any conditions of use.
7. Do not sit on or against or crush device—avoid sharp folds.
8. Never pull the device by the line cord. Do not use the cord as a handle.
9. The device does not include moisture resistant. Do not use while bathing or in the shower. Use only in a dry area. If any signs of deterioration (such as checking, blistering, or cracking) are visible, the device should not be used. Discard the device immediately.
10. Do not use on an infant or an animal.
11. Do not tamper with device in any way. There are no user-serviceable parts. If for any reason they do not function satisfactorily, return them to place of purchase. No repairs should be attempted by the consumer.
12. Do not use the device with liniment, salve, or ointment preparations that contain heat-producing ingredients. Skin burns could result.
13. Loop cord loosely when storing. Tight wrapping may damage cord and internal parts.
14. Do not put the device in the microwave.
15. Do not leave the device unattended, especially if children are present.
16. Do not use device underneath bedcover or other types of blankets.

17. Keep the device away from heated surfaces.
18. When not in use this device should be kept in a cool, dark, dry place, preferably in the
19. This is not a toy. Do NOT allow children to touch or play it.
20. Do not use the unit on an open wound, as heat may increase circulation and cause bleeding.
21. Do not use over thyroid gland (the neck area) if the user has been diagnosed with or has thyroid problems. No clinical data is available to support use for patients with thyroid gland
22. This product contains neoprene. If you are allergic to neoprene, do not use this product consulting a physician.
23. Use this appliance only for its intended use as described in this manual. Do not use at recommended by Conair.
24. Do not use outdoors.
25. Light may bring hurt to eyes. Do NOT look directly at the light of the LED.

## SAVE THESE INSTRUCTIONS!

### III. Warnings

1. The Conair® Flex Pain Relief is an electrical device. To avoid electrical shock:
  - a. Do not wash the unit with water.
  - b. Do not touch the plug or device with wet hands.
  - c. Do not touch unit to wet skin.
  - d. Do not use if cord or plug appears to be damaged.
  - e. Do not use the unit if it has been damaged, torn, or has come in contact with
  - f. Do not pull, lift, or carry unit by the power cord.
2. Do not operate the unit in areas where flammable or combustible products or vapors may be present. Explosion or fire may occur.
3. Do not allow children to use – this is not a toy. Do not leave small children unattended close to unit.
4. Not to be used by disabled, sleeping or unconscious persons.
5. Do not use with liniments, salves or ointments that contain heat-producing ingredients. A skin burn could result.

a cool, dark, dry place, preferably in the original carton. Do not use this product without first reading the instructions.

It may increase circulation and cause bleeding. If the user has been diagnosed with or suspects a thyroid gland problem, do not use this product without first reading the instructions.

Do not use attachments not described in this manual. Do not use attachments not described in this manual.

Use only with the supplied adapter. Do not use attachments not described in this manual.

# INSTRUCTIONS

Use only with the supplied adapter. Do not use attachments not described in this manual.

Do not use attachments not described in this manual.

Do not use attachments not described in this manual.

Do not use attachments not described in this manual.


Do not use attachments not described in this manual.

Do not use attachments not described in this manual.

## V. Product Details

The Conair® Flex Pain Relief comes completely assembled and ready to use. Your box will include the following:

1. Conair® Flex Pain Relief Flex Pad (soft neoprene)
2. Conair® Flex Pain Relief adjustable hook and loop support straps (2)
3. AC adapter with USB Port
4. User Guide and Operating Manual
5. ON/OFF Switch

## VI. Assembly Instructions

1. Attach hook and loop straps to back of pad.
2. Adjust hook and loop length to comfortably fit desired body area (back, leg, arm, etc.) for treatment.
3. Place Conair® Flex Pain Relief on body and secure hook and loop strap.
4. Insert USB plug from pad into any USB port or into the provided USB adaptor.
5. Plug adaptor into electrical outlet.
6. Press ON/OFF switch to the on position to start device.
7. Move ON/OFF switch to the off position to stop device.

## VII. Instructions for Use

1. Place Conair® Flex Pain Relief on body and secure hook and loop straps.
2. Insert USB plug from pad into any USB port or into the provided USB adaptor.
3. Plug adaptor into electrical outlet.
4. Press ON/OFF switch to the ON position to start device.
5. You can only see the Red LEDs; the others are in the infrared light spectrum and can't be seen with the human eye.
6. Treatment time is 20 minutes. The device may be used multiple times per day.
7. Move ON/OFF switch to the OFF position to stop device.

## VIII. Technical Specifications

- Model No: HP11FQ
- ONLY USE SUPPLIED USB adaptor.
- Device includes ON/OFF switch; all settings are set at the factory.
- The Conair® Flex is NOT user serviceable.
- Operating Temperature: 0°C to +27°C (32°F – 80°F)
- Relative Humidity: 85%
- Storage Temperature Range: -30°C to +55°C (-22°F to +130°F)
- Wavelengths: Infrared LED 880nm, RED LED 660nm
- Power: DC5V 1A

## IX. Cleaning Instructions

Clean LEDs monthly with mild soap and water. Do not use alcohol or other solvents.

## X. Frequently Asked Questions

- Q: How long before I can expect to see results?  
 A: Every person and every pain is different. Some people may see results immediately, while others may take several days to see results.
- Q: The system gets warm. Is it safe?  
 A: The infrared LEDs cause the heat. The system is designed to handle this.
- Q: Do the LEDs ever need to be replaced?  
 A: The LEDs are designed to last for the life of the device.
- Q: Are there other accessories available?  
 A: Yes. You can purchase additional accessories.

## XI. Disclaimer

This product is not intended to diagnose, cure, or prevent any disease. It is not a medical device.

## XII. Equipment

- Type BF Patient-Applicator
- Ordinary Equipment (not Class II)
- Equipment not suitable for use in the presence of flammable anesthetics
- Complies to EMC Standard EN 60601-1-2
- Complies to Medical Electrical Safety Standard IEC 60601-1

## IX. Cleaning Instructions

Clean LEDs monthly with rubbing alcohol and cotton cloth. Allow system to cool before applying rubbing alcohol. Allow alcohol to dry before use.

## X. Frequently Asked Questions

- Q: How long before I can expect results?  
 A: Every person and every pain is unique. Many people will experience reduced pain after their first use, for others it could take continual use for a week or more.
- Q: The system gets warm. Is that ok?  
 A: The infrared LEDs cause the system to warm the skin. This is normal, and the system is designed to handle this.
- Q: Do the LEDs ever need to be replaced?  
 A: The LEDs are designed to last years for the average consumer.
- Q: Are there other accessories for my system?  
 A: Yes. You can purchase accessories and replacement parts directly from Conair Corporation, Inc.

## XI. Disclaimer

This product is not intended to diagnose, cure or prevent disease. Conair Corporation, Inc. makes no claims, representations or warranties regarding the ability of this product to cure any physical, skin, or mental conditions. A qualified health professional should be consulted with regard to any condition requiring medical attention.

## XII. Equipment Classification

- Type BF Patient-Applied Part
- Ordinary Equipment (not protected against harmful ingress of moisture)
- Equipment not suitable for use with flammable anesthetics
- Complies to EMC Standard IEC 60601-1-2
- Complies to Medical Electrical Equipment Safety Standard 60601-1-1

## XIII. Warning Symbols

*Attention  
Consult  
Accompanying  
Documents*

The operator must read, understand, and follow all instructions in the accompanying documents, including all warnings, cautions, and precautions before using the medical device.

*Type BF  
Equipment*

This symbol indicates that the patient-applied part is Type BF (floating from electrical ground), which offers a specific level of safety.

*Standby/On*

This symbol indicates that the switch electronically cycles the DC power ON and OFF for part of the equipment.  
 NOTE: To disconnect the power, unplug from electrical outlet.

## XIV. Limited One-Year Warranty

Conair Corporation, Inc. warrants that your product (Conair® Flex Pain Relief) will be free from defects in materials and workmanship for one year from the date of purchase by the original owner and is not assignable or transferable to another party.  
 The product must be used in a manner as presented in this user guide and only with the provided power supply.

### ***What to Do if I Need This Warranty Service?***

To initiate warranty service, please contact Customer Service to receive an RMA (Return Merchandise Authorization) and further instructions. Proof of original purchase is required.

### ***How to Register Warranty?***

Conair, Inc. warrants that this product will be free from defects in materials and workmanship for 1 year from the date of purchase by the original owner. This warranty covers the original purchaser only and is not assignable or transferable to another party. Register your warranty online at [conair.com/registration](http://conair.com/registration)

## XV. FCC MARKINGS

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**Please register this product at  
[www.conair.com/registration](http://www.conair.com/registration)**

### **SERVICE CENTER CONAIR CORPORATION**

Service Department  
 150 Milford Road  
 East Windsor, NJ 08520

### **CONAIR CORPORATION**

Service Department  
 7475 N. Glen Harbor Blvd.  
 Glendale, AZ 85307

Visit us on the web at [www.conair.com](http://www.conair.com)

©2017 CONAIR CORPORATION  
 East Windsor, NJ 08520  
 Glendale, AZ 85307

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